

Product Datasheet

Name: TOX-IgM-p30

Description: Expressed in *E.coli* cell line. Blood sample from the tissue donors was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

Catalog No.	Source	Isotype	Unit	Buffer	Usage
P4103	<i>E.coli</i>	IgM	ml	0.01M PBS (pH7.4)	Detection

Format: Liquid

Concentration: ≥ 2 mg/ml

Purity: $\geq 90\%$ (SDS-PAGE)

Preservative: 0.09% sodium azide

Application: Each laboratory should determine an optimum working titer for use in its particular application.

Storage: Store at 2-8°C up to 21 days, -20°C for long term. Avoid multiple freeze-thaw cycles.

Note: Test systems and platform change may cause different results, please contact us if you have any questions.

For research or further manufacturing use only.

Nanjing Genstars Biotech Co.,Ltd.

Address: Layer 4, building A6&A7 Maple Science Park, Qixia District, Nanjing, China

Tel: +86-025-86793889

Email: marketing@genstars.com

Website: www.genstars.com



Website



Wechat