

## Product Datasheet

**Name:** HSV-2-IgG-gG II

**Description:** Expressed in *E.coli* cell line. Blood sample from the tissue donors was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

Catalog No.	Source	Isotype	Unit	Buffer	Usage
P2503	<i>E.coli</i>	IgG	mg	0.05M Tris (pH8.0)	Capture

**Analysis:** Recombinant Immunogen for antibody production

**Format:** Liquid

**Concentration:**  $\geq 2$  mg/ml

**Purity:**  $\geq 90\%$  (SDS-PAGE)

**Preservative:** 0.09% sodium azide

**Application:** Each laboratory should determine an optimum working titer for use in its particular application.

**Storage:** Store at 2-8°C up to 21 days, -20°C for long term. Avoid multiple freeze-thaw cycles.

**Note:** Test systems and platform change may cause different results, please contact us if you have any questions.

For research or further manufacturing use only.

### Nanjing Genstars Biotech Co.,Ltd.

Address: Layer 4, building A6&A7 Maple Science Park, Qixia District, Nanjing, China

Tel: +86-025-86793889

Email: [marketing@genstars.com](mailto:marketing@genstars.com)

Website: [www.genstars.com](http://www.genstars.com)



Website



Wechat